

Los Alamos National Laboratory is operated by the University of California for the United States Department of Energy under contract



S Alamos National Laboratory Los Alamos, New Mexico 87545



LA-9879-MS

UC-15

Issued: November 1984

The Assay of Plutonium Metal by Gamma Spectrometry and Calorimetry

David F. Bowersox Raymond P. Wagner





LOS Alamos National Laboratory Los Alamos, New Mexico 87545

Ъy

David F. Bowersox and Raymond P. Wagner

ABSTRACT

The nondestructive assay of plutonium buttons containing 1-3 kg plutonium agrees within 0.2% of the value calculated by chemical assay. These buttons were 4-15% 240 Pu by weight and contained 100 to 10000 ppm 241 Am. The nondestructive assay consisted of a gamma spectrometric determination of the isotopic distribution and a calorimetric determination of the power. Although further study is needed to extend the quantitative results to lower 240 Pu and 241 Am concentrations, the method is a rapid, nonintrusive technique for assaying plutonium.

I. INTRODUCTION

We are always interested in proposals for reliably assaying the plutonium content of materials in the Plutonium Facility at Los Alamos. If we could quickly determine an accurate and precise plutonium value, we could improve the present accounting procedures. The nondestructive assay of plutonium by gamma spectrometry and calorimetry is a promising method for such an assay. The isotopic distribution in typical material, containing 4-15% ²⁴⁰Pu, can be determined accurately in <4 h, and calorimetry can be completed in about the same time. Most shipping cans could be assayed as received with no need for unpackaging.

We determined the plutonium value for 400 metal buttons by gamma spectrometry/calorimetry, chemical assay/calorimetry, and chemical assay/mass measurement (plutonium-factor value). The ²⁴⁰Pu distribution value determined by gamma spectrometry was then used to select button blends for ingots of specified ²⁴⁰Pu composition. An independent value was that given by the shipper, which will be designated as the declared value. We do not know how the declared value was obtained; and, in fact, it may duplicate either the chemical/calorimeter or plutonium-factor value. We do not have a true plutonium value and do not claim that any one of the four values is more accurate than any other. Although differences in values for individual buttons are interesting, summarizing data is more useful. For

the purpose of this report, we will discuss and compare our determinations, the limitations of the present gamma instrument, and the accuracy and precision of the results.

II. EXPERIMENTAL

The gamma instrument is shown in Fig. 1. The intrinsic germanium detector, which is enclosed in lead shielding, is located directly under the sample chamber and beneath the glove box. The sample chamber consists of a 0.6-cm thick, 30-cm diam Lexan window in the glove box floor over the detector and a 2.5-cm-thick, 45-cm-diam lead disk with a variable diameter collimator (2.5-15.0 cm) directly over the window. Sample cans were placed on an adjustable platform over the lead disk. Both the diameter of the collimator and the sample-to-detector distance were varied to optimize counting statistics. The chamber was placed in a glove box to allow the assay of contaminated materials; however, we have used this system, up to now, only on packaged samples.

The electronics, analyzer, and computer for the instrument are approximately 10 m from the detector. A Canberra Series 80 Multichannel Analyzer (Fig. 2), interfaced with a PDP 11/34 computer and printer, is programmed to collect the data.

A sample that consisted of a double can containing the packaged button was placed on the platform in the sample chamber and the sample-to-detector distance adjusted to obtain 10000-12000 counts/s. At lower count rates, the can was placed on the lead shield (Fig. 3). The system was programmed to collect the peak areas at 25 gamma peaks between 125 and 415 keV for a preselected time, usually 3.5-12 h. Each peak is characteristic of the concentration of a given isotope, and by taking ratios of the areas of 238Pu peaks, 239Pu peaks, and 240Pu peaks to nearby 241Pu peaks and the ratio of the areas of 241Am peaks to 239Pu peaks, the distribution of the isotopes can be calculated. Since 242Pu has no known usable gamma peaks, its distribution was calculated from the other isotopes by isotopic correlation techniques. Fortunately, the ²⁴²Pu concentration is very small and its estimate can be quite poor without affecting the assay. The ratios, the fractions, the weight per cent of the plutonium isotopes, and the $^{241}\mathrm{Am}$ concentration were calculated by the computer. The results were printed and the spectra stored on a floppy disk. The gamma spectrometer was calibrated before these assays with a series of well-characterized plutonium oxides and a plutonium metal button. During the study, the metal button and a 1000-g can of well-characterized PuO2 were occasionally used to ensure that there were no changes in the calibration.




Fig. 1. Gamma detection station.




Fig. 2. Computer system.




Fig. 3. Can containing plutonium ready for gamma assay.




Fig. 4. Button with sample taken by drilling.

After unpacking, a sample was drilled out near the center of the bottom of each button and analyzed (Fig. 4). 'The ²³⁸Pu and ²⁴¹Am concentrations were determined by radiochemistry; the ²³⁹Pu, ²⁴⁰Pu, ²⁴¹Pu, and ²⁴²Pu distributions, by mass spectrometry. The methods were checked with standards from the NBS SRM series for isotopic analysis. The standard deviation in these analyses is <0.2%.² The sample for chemical assay must be representative of the entire button. A typical button is shown before cleaning in Fig. 5. Fifteen 5-g samples were taken from button HRA 147215 to check homogeneity.³ Twelve samples consisted of metal turnings from four main holes (Fig. 6). The other three samples consisted of pieces sheared off the button with a hydraulic chisel at many locations. The data are summarized in Table I. Clearly, the button is homogeneous, and the normal samples from the center of the bottom are representative.

The total plutonium in a button was calculated by dissolving portions of the drilled sample and assaying for total plutonium by controlled potential coulometry. The precision of this assay is >0.1%. By calculating grams of plutonium per grams of sample, multiplying by the button weight, and adding the plutonium in the residues, a plutonium value was obtained.

The procedure for determining the plutonium value was as follows:

- The packaged weight of the can and its contents were verified by reweighing.
- The plutonium isotopic distribution and the americium concentration were determined by gamma spectrometry.
- The total power of the packaged button was measured with a calorimeter.
- 4. The can was opened and the button was unpackaged, cleaned, and weighed. All loose oxide was collected and weighed. The plutonium in the packaging material was determined by neutron counting.
- 5. A small sample was drilled out near the center of the bottom of the button. The analytical group determined the isotopic distribution, the americium concentration, and the plutonium factor by assaying this material.

The specific power was determined from the isotopic distribution and americium concentration by a computer program.* The plutonium value was then

^{*}This information provided by T. E. Sampson of Los Alamos National Laboratory, December 29, 1981.




Fig. 5. Typical plutonium button before cleaning.




Fig. 6. Multiple sampling of button HRA 147215.

TABLE I
DATA FROM MULTIPLE SAMPLING OF HRA 147215

Sample ID		Weight Percent						Pu Value,g ^a	
						241 Am ,	Chem. Assay		
HRA 147+	238	239	240	241	242	bew	+Calor.	Pu Factor	
215	0.0081	94.09	5.68	0.202	0.021	1320	2140	2150	
006	0.0082	94.10	5.67	0.199	0.020	1330	2139	2147	
007	0.0081	94.08	5.69	0.202	0.020	1320	2139	2146	
010	0.0082	94.08	5.69	0.201	0.020	1370	2134	2147	
011	0.0085	94.09	5.68	0.200	0.020	1390	2131	2147	
012	0.0083	94.09	5.68	0.203	0.021	1370	2134	2148	
014	0.0083	94.08	5.69	0.201	0.020	1360	2135	2147	
015	0.0083	94.09	5.68	0.199	0.019	1360	2134	2148	
016.	0.0083	94.09	5.68	0.200	0.021	1350	2134	2136	
017	0.0084	94.10	5.67	0.200	0.020	1360	2135	2148	
018	0.0086	94.09	5.68	0.199	0.020	1390	2135	2143	
019	0.0079	94.09	5.68	0.203	0.021	1320	2141	2147	
020	0.0088	94.10	5.68	0.200	0.019	1330	2136	2147	
021	0.0081	94.09	5.68	0.200	0.020	1310	2141	2143	
022	0.0079	94.08	5.68	0.203	0.022	1300	2143	2150	
ř	0.0083	94.09	5.68	0.201	0.020	1350	2137	2146	
.s.d. (%)	±3.0	± 0.00%	±0.1%	±0.7%	±4.0%	± 2.1%	± 0.2%	± 0.2%	

^{*}SNM value is 2148 g.

calculated by dividing the total power by the specific power. This resulted in two independent values: one determined by calorimetry and gamma spectrometry; the other, by calorimetry and chemical analysis. A third plutonium value was obtained by multiplying the plutonium factor by the button weight and adding the plutonium in the packaging and that removed as oxide during button cleanup.

III. RESULTS

The approximate distributions of 240 Pu and 241 Am in each button are summarized in Table II. Plutonium values obtained by all three methods differed by <0.7%, as shown in Table III. The average value from three methods, plus the declared value for 409 buttons, was 798521 \pm 1396 g, which is equivalent to a standard deviation of 0.2%. The assay of 238 Pu in the first series that we examined, designated the EXCM, was higher by mass spectrometry due to contamination by 238 U. Radioassay determinations were used for 238 Pu on all subsequent chemical assays. The bias was apparently eliminated. The EXCM Series was not included in the totals because of this bias; however, the declared value and the values derived from calorimetry plus gamma spectrometry and from the plutonium factor are in good agreement.

The americium value determined by gamma spectrometry was consistently 10% higher than that determined by radioanalysis. This bias caused a low total plutonium value, which was corrected by dividing all gamma values by an instrument calibration factor of 0.982. We believe that this factor is due to the programming of this particular instrument and could probably be eliminated in another instrument. If the americium concentration was less than 400 ppm, one of the two peaks used to determine the ²⁴¹Am/²³⁹Pu ratio vanished. With our computer program, this caused an erroneously high americium result and, consequently, a low total plutonium value. We substituted a value of 400 ppm americium or used the radioassay value for americium in calculating the total plutonium in such cases, and because the instrumentation factor was largely due to high americium values, we did not apply the factor in such cases. The 240 Pu values at concentrations <3% were usually about 20% higher than normal with this instrument. We had not anticipated receiving material this low in 240Pu, and we did not have standards in this concentration range for calibration. In addition, for these buttons the 238 Pu concentrations were extremely low (0.005%) and the 241 Am concentrations also were <400 ppm. Our total plutonium values were low for these buttons; the calorimetry/chemical isotopic, plutonium-factor methods and declared values also

were not in good agreement. Therefore, the buttons containing <3% 240 Pu are not included in the summary of results. We believe that a gamma instrument programmed for this material would yield acceptable results.

TABLE II

THE APPROXIMATE DISTRIBUTIONS OF 240Pu AND 241Pu IN METAL BUTTONS

Series	Compo	sition
<u>ID</u>	240Pu (wt%)	²⁴¹ Am (ppm)
HRA 114	20-14	8000
HRA 99	10	5000
HRA 131	8-10	6000
ARF 147	6	1500
EXCM	6	1200
ARF 458	6	100-1500
ARF 469	6	<500
ARF 437	6	<500
ARF 484	6	<400
ARF 455	5-6	100-300
ARF 452	4-7	100-500
ARFSO	<2.5	<500
ARF 120A	<2.5	<500

TABLE III
SUMMARY OF RESULTS IN THE DETERMINATION OF Pu VALUES

Series	Number		Pu Value (<u>g)</u>			
<u>ID</u>	Buttons	Declared	Cal+Gamma	Cal+Chem	Pu Factor	Average	sd (g) (%)
HRA 114	27	46714	45410	47213	46700	46784	±301 (0.6)
HRA 99	25	88773	88870	88831	86816	88870	±112 (0.1)
HRA 131	26	53365	53393	53334	53236	53382	±151 (0.3)
ARF 147	60	124829	124313	124096	124576	124454	±318 (0.3)
EXCM ^a	21	68389	68614	67773	68498	68068	±869 (1.3)
ARF 458	57	108769	108042	108461	108785	108514	±348 (0.3)
ARF 469	35	67932	67116	67622	68235	67726	±478 (0.7)
ARF 437	49	94681	94439	93953	95151	94556	±449 (0.5)
ARF 484	12	24093	23773	23971	24134	23993	±162 (0.7)
ARF 455	23	43855	43503	43740	43997	43774	±209 (0.5)
ARF 452	23	44270	43811	44217	44477	44194	±279 (0.6)
HRA 12DA	52	102719	102844	102597	102597	102719	±140 (0.2)
Totals	409	800000	796814	798035	799235	798521	±1396 (0.2)

^aEXCM is not included in the totals because the ²³⁸Pu was determined by mass spectrometry, not radiochemistry. This caused low Pu values for the Cal+Chem method.

The average composition of the 45 buttons of the HRA 99 Series is 0.0310 wt% 238 Pu, 90.50 wt% 239 Pu, 8.79 wt% 240 Pu, 0.615 wt% 241 Pu, 0.75 wt% 242 Pu, and 4950 ppm americium. The isotopic distribution for the individual buttons, as determined by gamma spectrometry, is listed in Table IV. The precision of the assay increases with time and with the concentration of the isotope. For example, the estimated precision of the ²³⁸Pu measurement increases from 8.3% at 4.4 h to 2.5% at 14.0 h, and to 1.8% at 40 h for 0.03 wt% 238 Pu. The precision of the 240 Pu assay, based on a much weaker energy peak, increased from 4.2% to 2.1% to 1.2% over the same time period for buttons ~ 9 wt% 240 Pu. On the other hand, the precision of the measurements of both the 239 Pu and 241 Pu was >0.5% in <3 h. The average assay time for the HRA 99 Series was 10.3 h. The averages for the isotopic data, as determined from the gamma assay, are summarized in Table V. These include the isotopic distribution, the calculated power per gram of plutonium in the button, and the per cent of the specific power contributed by each isotope. The plutonium value is obtained by dividing the power determined by calorimetry by the specific power.

Because the per cent of the total power contributed by each isotope and the precision of each measurement are known, the precision of the contribution of each isotope to the calculation of the plutonium value can be estimated. The ²³⁸Pu contributes 5.61% to the total power even though the concentration is only 0.30%. The estimated precision of the ²³⁸Pu assay is equivalent to a precision of 0.30% in determining contribution to the total power. The major isotope in the buttons is ²³⁹Pu. Because the power generated by individual ²³⁹Pu atoms is far less per gram than that generated by either ²³⁸Pu or ²⁴¹Am, the contribution for the power, 55.75%, is much less than the concentration of the isotope. The estimated precision due to the ²³⁹Pu gamma assay is equivalent to 0.17% of the power.

The 240 Pu contributed 19.91% of the power and the estimated precision is 0.67%. The power from the 241 Pu is 0.67% and the precision is >0.01%. The 242 Pu does not contribute significantly to the total power. Finally, the contribution from 241 Am is 18.06% of the power with an estimated precision of 0.11%. If we assume that the precision of the calorimetric measurement is at least 0.3%, the estimated precision of the plutonium value by gamma spectrometry and calorimetry should be 0.8%. The precision of the ratio of the gamma/calorimetry derived plutonium value to the declared value is actually 0.7%.

TABLE IV

THE ISOTOPIC DISTRIBUTION OF THE HRA 99 SERIES

BY GAMMA SPECTROMETRY ASSAY

Los Alamos			D., T., 4.	N	n CA	
ID HRA 99+	238	239	240	pic Weight 241	242	
ILLA 991	230	239		241		ppm Am
JA 07	0.0267	90.94	8.42	0.5542	0.0620	4840
X 08	0.0264	91.56	7.83	0.5334	0.0548	4920
NX07	0.0189	90.09	9.26	0.5659	0.0710	5120
NX08	0.0241	90.46	8.86	0.5925	0.0706	5120
NX-9	0.0284	91.47	8.00	0.4584	0.0482	3720
			•			2,20
NXO9	0.0257	90.24	9.06	0.6052	0.0741	5260
NX11	0.0237	90.45	8.84	0.6152	0.0731	5420
NX12	0.0296	91.06	8.25	0.6032	0.0660	5080
NX13	0.0261	90.47	8.83	0.5993	0.0711	5190
NX14	0.0283	90.85	8.46	0.5989	0.0675	5070
NX15	0.0408	89.12	9.93	0.7938	0.1092	6 520
NX16	0.0459	89.04	9.91	0.8842	0.1215	7160
NX17	0.0464	88.55	10.35	0.9138	0.1327	7410
NX18	0.0471	88.44	10.51	0.8801	0.2301	7240
NX19	0.507	88.97	9.94	0.9144	0.1264	7420
NX20	0.0504	88.74	10.17	0.9138	0.1298	7310
NX21	0.0481	88.35	10.56	0.9139	0.1360	7420
NX22	0.0511	88.32	10.55	0.9540	0.2060	7700
NX24	0.0250	91.38	7.99	0.5456	0.0547	3610
NX35	0.0267	91.03	8.35	0.5389	0.0597	4190
NX36	0.0276	90.90	8.48	0.5362	0.0605	4160
NX 37	0.0261	91.36	8.02	0.5377	0.0569	4280
NX38	0.0258	90.53	8.86	0.5295	0.0629	4190
NX39	0.0341	92.04	7.41	0.4742	0.0456	4410
NX40	0.0266	90.43	8.95	0.5350	0.0644	4250
NX41	0.0271	90.82	8.53	0.5572	0.0634	4370
NX42	0.0285	90.32	9.02	0.5643	0.0687	4450
NX43	0.0328	90.11	9.10	0.6791	0.0837	5360
NX44	0.0289	90.30	8.99	0.6084	0.0738	4700
NX45	0.0300	90.36	8.97	0.5696	0.0688	4530
357700	0.0070	01.06	2 55	0.5050	0.0400	2000
MX80	0.0273	91.86	7.55	0.5053	0.0498	3000
MX81	0.0245	91.51	7.83	0.5812	0.0598	6470
NX86	0.0254	90.67	8.71	0.5364	0.0625	4150
NX87	0.0273	91.13	8.26	0.5305	0.0580	4060
NX88	0.0260	90.79	8.59	0.5314	0.0609	4050
NX89	0.0276	90.84	8.54	0.5309	0.0605	3940
NX90	0.0278	90.76	8.62	0.5305	0.0611	4090
NX91	0.0257	90.70	8.48	0.5295	0.0597	4120
NX92	0.0272	91.00	8.39	0.5279	0.0588	4000
NX93	0.0253	90.65	8.73	0.5367	0.0627	4110
en))	C.32.0	30.03	0.75	0.5501	0.0021	7110
XX94	0.0266	90.30	9.03	0.5727	0.0698	4350
MX95	0.0287	89.65	9.65	0.5928	0.0782	4650
W X98	0.0266	92.50	6.96	0.4685	0.0419	2910
X 01	0.0461	91.96	7.34	0.6028	0.0575	4030
MX99	0.0275	91.06	8.31	0.5353	0.0590	4270
				+.000	3.0370	72,0

TABLE V
ISOTOPIC DATA AVERAGES FOR THE HRA 99 SERIES

			Specific	
[sotope	Isotope (W/g)	Isotopic Distribution	Power (W/g) Pu	Per Cent of Total Power
238	0.5672	0.03098%	1.7569×10^{-4}	5.61
239	0.001929	90.50%	17.46 × 10 ⁻⁴	55.75
240	0.007098	8.786%	6.236×10^{-4}	19.91
241	0.003390	0.6152%	0.2085 x 10 ⁻⁴	0.67
242	0.0001159	0.07584%	0.0009×10^{-4}	0.00
Am-241	0.1142	4950 ppm	5.654×10^{-4}	18.06
Total			31.3253 x 10 ⁻⁴	100.00

TABLE VI

THE RATIOS OF ISOTOPIC VALUES; GS/CHEM FOR THE HRA 99 SERIES

Los Alamos	Ratio: Gamma Assay/Chemical Assay							
ID K RA 99+	238	239	240	241	242	Am		
JA07	0.988	0.995	1.059	1.008	0.775	1.093		
X08	0.912	0.999	1.022	0.994	0.481	1.076		
X07	1.129	1.002	0.984	1.006	1.044	1.008		
1X08	1.053	1.002	0.986	1.016	1.008	1.057		
IX-9								
1.4-9	1.208	0.997	1.039	0.987	1.461	1.191		
IX09	0.991	0.997	1.035	1.015	1.029	1.185		
TX 1 1	0.912	1.001	0.991	1.007	0.897	1.058		
IX12	1.095	1.004	0.959	1.010	0.970	1.047		
TX 13	0.966	0.997	1.029	1.011	1.046	1.053		
X14	1.053	1.002	0.978	0.999	0.937	1.074		
X15	1.020	0.998	1.022	1.009	0.628	1.042		
X16	0.997	1.004	0.971	1.013	0.552	1.065		
X17	0.967	1.001	0.997	1.009	0.561			
X18						1.083		
	1.025	0.998	1.028	1.008	0.593	1.081		
X19	1.056	1.006	0.962	1.018	0.543	1.085		
X20	1.050	1.003	0.984	1.013	0.557	1.069		
X21	1.023	0.999	1.017	1.004	0.575	1.076		
X22	1.022	1.002	0.993	1.015	0.792	1.059		
X24	0.967	1.002	0.980	1008	0.545	1.074		
X35	1.068	1.001	0.987	1.002	1.328	1.096		
X 36	1.028	1.000	1.002	1.002	1.345	1.066		
X37	1.019	1.005	0.949	1.000	1.264	1.079		
X38	0.968	0.995	1.048	0.993	1.430	1.115		
X39	1.004	1.000	1.004	1.000				
X40	1.030	0.995	1.054	0.992	0.449 1.431	1.103 1.102		
A-10	1.030	0.333	1.054	0.332	1.431	1.102		
X41	1.062	1.000	0.995	1.000	1.321	1.081		
X42	1.032	0.996	1.045	0.994	1.164	1.109		
X43	1.066	1.002	0.983	1.005	1.087	1.078		
X 44	1.006	O.998	1.017	0.998	1.229	1.087		
X 45	1.077	0.997	1.033	0.979	1.324	1.032		
x8 0	1.008	0.996	1.058	0.995	0.754	1.076		
X81	0.720	1.000	1.003	0.993	1.107	1.102		
X86	0.999	0.997	1.034	1.005	1.421	1.084		
X87	1.092	1.001	0.987	1.006	1.350	1.069		
X88	1.040	0.998	1.026	1.010	1.128	1.003		
von	1 000	0 000	1 000	1 000	1 /07			
X89	1.098	0.899	1.020	1.009	1.407	1.042		
X90	0.995	0.998	1.024	1.002	1.421	1.100		
X91	1.046	0.999	1.013	1.010	1.389	1.132		
X92	1.086	1.000	1.002	1.001	1.368	1.056		
X93	0.998	0.997	1.033	1.005	1.394	1.036		
X94	0.994	0.997	1.033	1.005	1.316	1.054		
X95	1.026	0.992	1.079	1.002	1.397	1.058		
X98	0.989	1.000	1.006	1.000	0.607	1.133		
X0 1	1.020	0.998	1.040	1.002	0.492	1.064		
X99	1.069	1.001	0.984	0.998	1.312	1.044		
ī	1.021	0.999	1.011	1.004	1.027	1.079		
.d.(%)	±0.07	±0.003	±0.029	±0.008	±0.345	±0.035		

The ratios of the isotopic values determined by gamma spectrometry to those determined by chemical assay are given in Table VI. The ²³⁸Pu by gamma assay averages 2.1% high, which is equivalent to a 0.1% lower total plutonium value. The ²³⁹Pu ratio is 0.999 so that the average increase in the plutonium value by gamma assay is <0.1% (0.056%). The ²⁴⁰Pu by gamma assay is 1.1% higher than by mass spectrometry, which represents a 0.2% lower total plutonium value. The ²⁴¹Pu and ²⁴²Pu ratios (<0.01%) are not large enough to cause significant changes in the plutonium value. The major differences are due to the high bias in the americium value, 7.9%. This will cause the average plutonium value calculated by the gamma data to be 1.4% lower than that derived from the chemical assay. The total calculated difference in the plutonium value is 1.71% lower for the gamma assay for this series of buttons. In determinations for a large number of buttons, the average was, indeed, 1.82% lower than the declared value. The difference between calculated and actual bias may be due to a small additional bias in the calorimetry.

The plutonium values for the HRA 99 Series are summarized in Table VII. Although individual assays vary, the total plutonium values are in good agreement. These totals are also within 0.2% of the declared value. Where the gamma/calorimetry values have been corrected for instrument bias (1.018 times the raw value), the three methods are equivalent.

IV. CONCLUSIONS

The calorimetry/gamma spectrometry method for determining the plutonium content is attractive because it is rapid, nondestructive, and nonintrusive. If the instruments are calibrated with similar standards, if the same program is employed, and if reported data are standardardized, this technique will provide an excellent procedure for shipper-receiver determination of plutonium value. At Los Alamos, Group Q-1 has completed a second gamma system incorporating the changes indicated by this study. It can be used for routine checks of plutonium content in buttons, oxide, and scrap. If it is used as planned, it should improve resolution of the plutonium content of all plutonium-containing materials.

ACKNOWLEDGMENTS

S. T. Huse, T. E. Sampson, J. L. Parker, and S. S. Johnson prepared and calibrated the gamma spectrometer system for this study and also gave helpful advice. William Haag, Cynthia Heinberg, Harriet Muffly, Phyllis Newlander, and Jimmie Torres performed many of the gamma assays, calculated the specific power, and


TABLE VII
SUMMARY OF Pu VALUES FOR THE HRA 99 SERIES

Los Alemos	İ	Pu Value,		Ratios			
IL	l	Calor.+	Calor.+	Pu Factor +	Games	Ches	Pu Facto
HRA 99+	Declared	Gamma	Ches.	Oxide + n.c.	Declared	Declared	Declared
JAD7	920	921	927	922	1.001	1.008	1.002
MOR	1775	1786	1775	1771	1.006	1.000	0.998
WX07	2010	2021	1997	2010	1.006	0.994	1.000
MXO8	1945	1959	1952	1937	1.007	1.004	0.996
3/X-9	2059	2027	2071	2057	0.985	1.006	0.999
#X09	1852	1853	1892	1853	1.000	1.022	1.000
10X11	2006	2019	2000	2006	1.006	0.997	1.000
NX12	2041	2053	2033	2041	1.006	0.996	1.000
NX13	2020	2017	2012	2019	0.999	0.996	1.000
MX14	1927	1948	1944	1921	1.011	1.009	0.997
MX15	2035	2026	2026	2033	0.996	0.996	0.999
MX16	2000	2004	1997	1994	1.002	0.999	0.997
NX.17	1979	1994	199	1976	1.008	1.007	0.999
MX1B	1977	1992	1923	1978	1.008	0.973	1.001
MX19	2004	2011	2020	2003	1.004	1.008	1.000
MX20	2047	2050	2056	2032	1.002	1.004	0.993
NX21	2031	1986	2042	2030	0.978	1.004	1.000
#D.22	1992	1987	1989	1974	0.998	0.999	0.991
NX24	2032	2053	2029	2020	1.010	0.999	0.991
MX35	1930	1951	1949	1934	1.011	1.010	1.002
ددسن	1,750	1731	1747	****	1	1.010	1.002
#X36	1941	1943	1936	1938	1.001	0.997	0.999
#X37	1998	2014	1991	1997	1.008	0.997	1.000
90 .38	1917	1910	1920	1918	0.996	1.002	1.001
WX39	1953	1949	1950	1955	0.998	0.999	1.001
EX-0	1778	1769	1782	1776	0.995	1.002	0.999
#X41	1949	1956	19 51	1946	1.004	1.001	0.999
ID:42	1907	1896	1910	1903	0.994	1.002	0.998
IIX43	1841	1844	1844	1838	1.002	1.002	0.998
MX44	1969	1970	1969	1969	1.000	1.000	1.000
MX45	2092	2068	2072	2092	0.989	0.990	1.000
WX80	2081	2072	2070	2083	0.996	0.995	1.001
MX81	269	270	267	269	1.004	0.993	1.000
PIX86	2063	2070	2075	2057	1.003	1.006	0.997
XX67	2070	2073	2066	2069	1.003	0.998	1.000
MX88	2045	2053	2055	2039	1.001	1.005	
MADO	2043	2033	2055	2039	1.004	1.005	0.997
MX89	1993	2003	2001	1977	1.005	1.004	0.992
90X90	2049	2058	2063	2044	1.004	1.007	0.998
30 00 1	2029	2036	2050	2023	1.003	1.010	0.997
EX92	2066	2081	2074	2066	1.007	1.004	1.000
MX93	2052	20 52	2041	2042	1.000	0.995	0.995
XX 94	2063	2065	2059	2057	1.001	0.998	0.997
XX95	2059	2032	2047	2057	0.987	0.994	0.999
MX98	2063	2057	2051	2063	0.997	0.994	1.000
XO1	2038	2038	2017	2036	1.000	0.990	0.999
MX99	1903	1903	1885	1903	1.000	0.991	1.000
_	1				l ——		
Σ	86770	86840	86828	86613	1.001	1.000	0.999
AX		70	58	-157	s.d. ±.007	±.008	±.003
AX,%		0.1	0.1	- 0.2	s.d. 07%	0.8%	0.3%

X₆= 86763 a.d.(%) ± 104 (0.1) tabulated much of the data. The metal fabrication section of Group CMB-11 unpackaged the ingots and sampled them; the count room of Group CMB-11 provided the calorimetry; and the analytical group, Group CMB-1, did all the conventional assays. We wish to thank them all for their help and encouragement.

REFERENCES

- 1. S. T. Hsue, T. E. Sampson, J. L. Parker, S. S. Johnson, and D. F. Bowersox, "Plutonium Isotopic Composition by Gamma-Ray Spectroscopy," Los Alamos Scientific Laboratory report LA-8603-MS (November 1980).
- 2. RDT F11-1 "Analytical Chemistry Methods for Mixed Oxide Fuel," Method 5, Plutonium and Uranium Isotopic Composition, 1972.
- 3. R. P. Wagner and J. D. Torres, "Multiple Sampling of a Metal Button," CMB-11-NMO-82-001, January 4, 1982.
- 4. G. R. Waterbury et al., "Controlled-Potential Coulometric and Potentiometric Titrations of Uranium and Plutonium in Ceramic-Type Materials," Los Alamos Scientific Laboratory report LA-4537 (October 1970).



Los Alamos